

# Injection Trainer

Part No: 00310

**Learning Has Never  
Been More Lifelike**

Epidermis layer may be peeled back, to release the subcuticular injection liquid, and replaced in seconds

Multiple tissue layers: epidermis, dermis, fat & muscle

Trainer is easy to assemble and disassemble

All components are latex free

Trainer can be strapped to a trainee's arm or thigh for role playing


Liquid injected into the pad can be squeezed out



# Injection Trainer

Part No: 00310

Multi-layered soft tissue pad for the practice of intradermal, subcutaneous, and intramuscular tissue injection techniques.



## Skills Gained:

- Professional-to-patient communication
- Subcutaneous injection
- Intradermal injection and Intramuscular injection
- Management of tissue

## Product Qualities:



- Multiple tissue layers: epidermis, dermis, fat & muscle
- Long-lasting replaceable epidermis
- Trainer can be strapped to a trainee's arm or thigh for role playing



- Epidermis layer may be peeled back, to release the subcuticular injection liquid, and replaced in seconds
- Liquid injected into the intramuscular pad can be squeezed out
- Trainer is easy to assemble and disassemble



- All components are latex free



## You'll Receive:

- 00311 1 Skin Pad & Muscle Block
- 00320 1 Epidermis for Injection Trainer (Pack of 2)
- 1 Strap-on Base Plate
- 1 Soft Tissue Tray with Liner

Packed weight (kg) & dimensions (mm) 600 g 210 l 160 w 110 h



Learn more at

[limbsandthings.com](https://limbsandthings.com)

Limbs & Things Ltd. Sussex Street, St Philips, Bristol, BS2 0RA, UK

T: 00 44 (0)117 311 0500 F: 00 44 (0)117 311 0501

E: UK Customers - [gbsales@limbsandthings.com](mailto:gbsales@limbsandthings.com)

E: International Customers - [sales@limbsandthings.com](mailto:sales@limbsandthings.com)